



IV Administration Set With GVS Speedflow[®] 0.2 μ m Filter

Our IV Administration Set with a 0.2 μ m filter consists of a universal spike, male swivel luer lock, roller clamp, drip chamber with a 15 μ m filter, and a 0.2 μ m in line filter. Each set is individually blister-packaged in quantities of 50 per box.

- 20 drops/mL
- One (1) Y-injection site 6" from distal end
- Male swivel luer lock connector
- Universal vented/non-vented spike
- Roller clamp
- 15 μ m filter located in drip chamber
- GVS Speedflow[®] 0.2 μ m Filter
- 92" length
- DEHP-Free
- Latex-Free



Product # TCBINF022AE

Visit www.tcbiomedix.com and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.