

6356 Manor Lane
South Miami, FL, 33143
Phone 866.593.8444
Fax 305.590.5537
info@tcbiomedix.com

TrueCare^{USA}
BIOMEDIX



IV Administration Set with Flow Regulator

- 100" length
- 2-Needle Free Y- Sites
- Universal Vented/Non-Vented Spike
- 10-drops/ml with a 15-micron filter in the drip chamber
- Roller Clamp
- Male Swivel Luer

Backcheck valve

Flow Regulator that delivers up to 300 ml per hour

DEHP Free PVC and Latex Free

Priming volume ~19ml

40 sets/Box



Visit www.tcbiomedix.com and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.

Product # TCBINF2NFRB
MKT-022 Rev