



## Pharmacy Transfer Tubing with Large Bore MLL

*Our custom Single Lead Pharmacy Fluid Transfer Set features a universal male luer lock and silicone tubing insert.*

- For use with Baxa Repeater® Pump
- Large bore intake tubing
- Universal male luer lock
- Silicone tubing insert
- DEHP-Free
- Latex-Free
- Pyrogen-Free
- Sterile



Product # TCBPTT021

Visit [www.tcbiomedix.com](http://www.tcbiomedix.com) and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.