



## PharmAssist<sup>®</sup> Compatible Solution Transfer Tubing

*Our custom PharmAssist<sup>®</sup> compatible transfer tubing pump sets are designed for use with the Excelsior PharmAssist<sup>®</sup> pump and are constructed for high volume use.*

- Universal vented spike
- Silicone tubing
- DEHP-Free
- Latex-Free
- One pinch clamp
- Male luer lock at distal end
- High volume capacity
- Pyrogen-Free



Visit [www.tcbiomedix.com](http://www.tcbiomedix.com) and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.

Product # TCBPTPS1