

6356 Manor Lane
South Miami, FL, 33143
Phone 866.593.8444
Fax 305.590.5537
info@tcbiomedix.com



Sterile Lab Gowns

TrueCare offers two different types of disposable sterile cleanroom gowns. Our sterile SMS 40 gram weight white, anti-static treated cleanroom gown has a stand-up collar with a Velcro closure in the back and elastic wrist cuffs. These gowns are individually packaged and come 50 to a carton. Our sizes consist of small, medium and large.

Our Sterile Chemo Compounding Gowns feature Barrier Max, ISO-certified breathable fabric. The gowns are impervious, low linting and Latex-free. Our gowns offer superior protection and comfort in the safe handling of hazardous, chemotherapy and other cytotoxic drugs. These gowns are individually packaged and come 50 to a carton. Our sizes consist of small, medium and large.

Product # TCBA54ST-S
Product # TCBA54ST-R
Product # TCBA54ST-XL

- Sterile Clean Room Gown, 52gsm+/-3gsm Breathable PE/PP, Level 3 impervious, CHEMO TESTED stand up collar with velcro closure in the back, with elastic wrists, 1 pc/sterile bag, 50 pcs/ctn

Product # TCBA40ST-S
Product # TCBA40ST-R
Product # TCBA40ST-XL

- Sterile Clean Room Gown, 40gsm+/-3gsm White SMS, anti-static treated, stand up collar with velcro closure in the back, with elastic wrists, 1 pc/sterile bag, 50 pcs/ctn



Visit www.tcbiomedix.com and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.

Product # TCBA40ST-S
Product # TCBA40ST-R
Product # TCBA40ST-XL
Product # TCBA54ST-S
Product # TCBA54ST-R
Product # TCBA54ST-XL