



Sterile Cleanroom Wipers

Our sterile cleanroom wipers are Sterile Gamma Irradiated to a 10⁻⁶ Sterility Assurance level. They are constructed of an extremely low linting 68gm./sq. meter 55/45 cellulose blend. They are appropriate for critical areas up to and including ISO5 Cleanroom environment. The Cleanroom Wipers feature maximum absorbencies of 20-25 ml., and are available in 9" x 9" and 12" x 12" sizes, and available in different quantities.

Product # TCBWIP09

- 300 wipers per pack, 12 packs per case, a total of 3600 wipers per case

Product # TCBWIP09SP-20

- 20 wipers per bag, 15 bags per bundle, 12 bundles per case, a total of 3600 wipers case

Product # TCBWIP12

- 150 wipers per pack, 8 packs per case, a total of 1200 wipers per case

Product # TCBWIP12SP-20

- 20 wipers per bag, 8 bags per bundle, 8 bundles per case, a total of 1280 wipers per case

