



## Micro Bore IV Extension Set

*Our Micro Bore IV Extension Sets are equipped with male and female luer locks and a slide clamp. The sets are DEHP-Free and Latex-Free, and are packaged in quantities of 50 per box.*

### Product # TCBEXT001

Male and female luer locks, slide clamp, 7" length, DEHP-Free, Latex-Free, 0.21mL priming volume

### Product # TCBEXT003

Male and female luer locks, slide clamp, 15" length, DEHP-Free, Latex-Free, 0.35mL priming volume

### Product # TCBEXT004

Male and female luer locks, slide clamp, 30" length, DEHP-Free, Latex-Free, 0.6mL priming volume



Product # TCBEXT001  
Product # TCBEXT003  
Product # TCBEXT004

Visit [www.tcbiomedix.com](http://www.tcbiomedix.com) and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.