

6356 Manor Lane
South Miami, FL, 33143
Phone 866.593.8444
Fax 305.590.5537
info@tcbiomedix.com



IV Administration Set

With 0.2 GVS SpeedFlow® Micron Filter
and GVS EasyDrop® Flow Regulator

Our IV Administration Set with 0.2 micron filter and flow regulator consists of a universal spike, male swivel luer lock, roller clamp, drip chamber with 15 micron filter, and 0.2 micron in-line filter. Each set is individually blister-packaged in quantities of 40 per box and features patient information directly on the exterior of the package.

- 20 drops/ml
- One (1) Y-injection site 6" from distal end
- Male swivel luer lock connector
- Universal vented/non-vented spike
- Roller clamp
- 15 micron filter located in drop chamber
- 0.2 micron GVS Speedflow® in-line filter
- GVS EasyDrop® Flow Regulator
- 92" length
- DEHP-free
- Latex-free



Product # TCBINF044G