



# TCBAMKST

## Sterile Face Mask

- + Excellent bacterial and particle filtration efficiency of  $\geq 98\%$ .
- + Fully enclosed, soft bendable nosepiece with knitted elastic ear loops.
- + External and internal layer of polypropylene spunbond, middle layer of polypropylene melt blown.
- + Recommended for ISO Class 5 and higher cleanroom environments.



Product #	Size	Features	Quantity
TCBAMKST	6.8" x 3.5"	<ul style="list-style-type: none"> <li>• White 3-ply mask with ear loops</li> <li>• Comfortable and Breathable</li> <li>• Meets ASTM F21000-20 Level 3 Requirements</li> <li>• BFE and PFE <math>\geq 98\%</math></li> <li>• Differential Pressure &lt; 6.0</li> <li>• Splash Resistance Pressure &gt; 160 mmHg</li> <li>• Flammability 16 CFR 1610: Class 1</li> <li>• Latex-Free</li> <li>• EO sterilized to <math>10^{-6}</math> Sterile Assurance Level</li> </ul>	Individually packed 20 pcs/bag 20 bags/case 400 pcs/case

