



## Sterile Lab Gowns for Hazardous Drug Compounding

Our Sterile Cleanroom Gown is disposable and composed of 52gsm weight anti-static material designed to meet USP 800 requirements. The Gown is breathable with PE/PP and offers Level 3 protection against chemo drugs. They come individually packaged, 1pc/sterile bag, 50 bags/case.

- Polypropylene spunbond fabric with a polyethylene exterior coating
- ASTM tested
- ISO-certified breathable fabric
- Level 3 protection
- Low linting
- Thumb loops
- Seam seals
- Stand up collar
- Elastic cuffs
- Latex-Free
- Available in sizes Small, Regular, XL and XXXL



Visit [www.tcbiomedix.com](http://www.tcbiomedix.com) and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets.

Product # TCBA54ST-S  
Product # TCBA54ST-R  
Product # TCBA54ST-XL  
Product # TCBA54ST-XXXL