




TCBEXT0.2

Filtered IV Extension Set 0.2µm

- + Engineered to provide advanced filtration, eliminating endotoxins, bacteria, and air bubbles from infusions.
- + Compact extension set is ideal for applications requiring precise filtration, offering exceptional reliability and protection in various clinical environments.
- + Minimizes the risk of infection and contamination, making it ideal for medications requiring filtration.
- + Reduces tangling and line adjustments, ensuring a smooth infusion experience in small or limited spaces.



Product #	Length	Priming Vol.	Features	Quantity
TCBEXT0.2	12"	4.5mL (approx.)	<ul style="list-style-type: none"> • PES Filter Membrane - 0.2µm • Male and Female Luer • Slide Clamp • Standard Bore • DEHP-Free • Latex-Free • Pyrogen-Free 	50 per box

