




# TCBEXT001

## 7" Mini-Bore Extension Set

- + Allows for easier clinician use.
- + Can be saline locked.
- + Reduce micro-movements at the IV insertion sites.
- + When using the clamp immediately prevent the backflow of blood, reducing exposure of nurses and other healthcare staff.



Product #	Length	Priming Vol.	Features	Quantity
TCBEXT001	7"	0.21mL (approx.)	<ul style="list-style-type: none"> <li>• Male and Female Luer</li> <li>• Slide Clamp</li> <li>• DEHP-Free</li> <li>• Latex-Free</li> <li>• Pyrogen-Free</li> </ul>	50 per box

