




TCBEXT1.2

Filtered IV Extension Set 1.2µm

- + Designed to enhance safety and efficacy in IV therapy by removing particles, bacteria, and other contaminants before they reach the patient.
- + Ideal for infusing lipids, nutrients, and medications.
- + 12" length is optimized for clinical convenience, allowing for easier management and reduced tubing bulk, making it ideal for tight spaces or precise line extensions.
- + Enables seamless integration into existing IV setups without additional adapters, saving time and ensuring a dependable connection.



Product #	Length	Priming Vol.	Features	Quantity
TCBEXT1.2	12"	4.5mL (approx.)	<ul style="list-style-type: none"> • PES Filter Membrane - 1.2µm • Male and Female Luer • Slide Clamp • Standard Bore • DEHP-Free • Latex-Free • Pyrogen-Free 	50 per box

