




TCBEXTNFBS-7

7" Standard Bore Bonded Needle Free Extension Set

- + Used for enhanced safety, efficiency, and reliability.
- + Lower chance of accidental disconnect.
- + Eliminates a connection point, reducing the potential for bacterial entry and contamination.
- + Eliminates the need for manual assembly, saving time in critical care settings.
- + Standard bore tubing for higher flow rates makes this beneficial in fluid resuscitation, blood transfusions, and high volume medication delivery, where rapid administration is critical.
- + Minimizes resistance, ensuring consistent and reliable flow, particularly when used with gravity fed IV systems.



Product #	Length	Priming Vol.	Features	Quantity
TCBEXTNFBS-7	7"	1.5mL (approx.)	<ul style="list-style-type: none"> • Slide Clamp • DEHP-Free • Latex-Free • 7 Inch Standard Bore • Bonded Needle-Free Valve • Rotating Male Luer Lock • Sterile 	50 per box

