




TCBEXTNFBM-14

14" Mini-Bore Bonded Needle Free Extension Set

- + Used for enhanced safety, efficiency, and reliability.
- + Lower chance of accidental disconnects or leakage.
- + Eliminates a connection point, reducing the potential for bacterial entry and contamination.
- + Eliminates the need for manual assembly, saving time in critical care settings.
- + Use for precision drug deliveries.
- + Longer extension allows patients to move more freely without pulling on the IV site.
- + Helps avoid frequent catheter manipulation, which is beneficial for oncology, ICU, and long-term therapy patients.



Product #	Length	Priming Vol.	Features	Quantity
TCBEXTNFBM-14	14"	0.6mL (approx.)	<ul style="list-style-type: none"> • Slide Clamp • DEHP-Free • Latex-Free • 14 Inch Mini-Bore • Bonded Needle-Free Valve • Rotating Male Luer Lock • Sterile 	50 per box

