



TCBWALC30

Sterile Isopropyl Alcohol Wipers

- + Low linting-ideal for use in pharmaceutical, biotechnology, and electronics industries.
- + Double bagged for added cleanroom protection.
- + Spunlace nonwoven (60 viscose - 40 PET).
- + Provides ample coverage for efficient cleaning and disinfection.
- + Presaturated with (70% isopropyl alcohol) for efficient removal of contaminants and residues.
- + Evaporates quickly, leaves no residue.
- + Resealable packaging - maintains moisture retention for long-lasting effectiveness.



Product #	Size	Features	Quantity
TCBWALC30	9" x 9"	• Gamma Irradiated to a 10⁻⁶ Sterility Assurance Level • Saturated with 70% USP / 30% DI Water • Alcohol is Filtered Twice Through 0.2µm Filters and then Terminally Sterilized • Extremely Low Linting, 65gsm Cellulose Blend • Double Bagged for Added Cleanroom Protection • Isopropyl Alcohol	30 wipers per pack 24 packs per case

