



TCBA54ST

Sterile Lab Gowns for Hazardous Drug Compounding

- + Composed of 52gsm $\pm$ 3gsm weight anti-static material.
- + Breathable fabric with tie-back, closed front and velcro collar.
- + Tape reinforced impervious seams.
- + Meet USP 800 requirements.
- + Hazardous drug compatible.
- + ASTM F3267-22 testing data available upon request.



Product #	Size	Features	Quantity
TCBA54ST-S	Small	• Polypropylene Spunbond Fabric with a Polyethylene Exterior Coating • ASTM F3267-22 Tested • Breathable Fabric • Level 3 protection per AAMI PB 70 • Low Linting • Thumb Loops • Seam Seals • Stand Up Collar • Elastic Cuffs • Latex-Free • Sterile to 10⁻⁶ Sterility Assurance Level	1pc per sterile pouch 50 pouches per case
TCBA54ST-R	Regular		
TCBA54ST-XL	XL		
TCBA54ST-XXXL	XXXL		

