



TCBCOMBICAP

Red Caps

- + Safe and secure connections to luer lock syringes.
- + Made from durable high-density polyethylene.
- + Offers security and versatility when capping syringes.
- + Works as catheter cut-off.
- + Designed to securely fit standard IV tubing and syringes.
- + Packaged individually to maintain sterility and prevent cross-contamination.



Product #	Features	Qty/Unit
TCBCOMBICAP	• A Male Luer Lock on One End and a Female Luer Lock on the Other • Individually Packaged • Sterile	100 per box

