



TCBPTT001

Pharmacy Transfer Tubing

- + Ensures seamless compatibility with various IV fluid containers, minimizing the need for multiple set types and streamlining pharmacy workflow.
- + Provides precise, consistent fluid delivery, optimizing accuracy in high-volume sterile compounding.
- + Allows efficient filling of multiple IV bags or syringes from a single bulk container, reducing manual handling time and increasing productivity in high-volume settings.



Product #	Features	Quantity
TCBPTT001	• Compatible with Baxa Repeater® Pump • Single Lead • Universal Vented Spike • Silicone Tubing Insert • Male Swivel Luer Lock at Distal End • DEHP-Free, Latex-Free, and Pyrogen-Free • Large Bore Tubing	10 per box

