

TCBPTT062

Bifurcated Pharmacy Transfer Tubing

- + Optimized for high-throughput pharmacy compounding supports efficient sterile compounding by minimizing bottle hangeovers and streamlining workflow.
- + Bifurcated large volume intake tubing enables simultaneous transfer of fluids from two large-volume sources, optimizing efficiency and reducing handling time.
- + Specifically designed with silicone tubing insert for seamless compatibility with the Baxa Repeater[®] Pump, fluid transfer with consistent flow rates and minimal disruptions.



Product #	Features	Quantity
TCBPTT062	• For use with Baxa Repeater[®] Pump • Two Large Volume Intake Tubing • Universal Vented Spike • Silicone Tubing Insert • Male Luer Lock at Distal End • Sterile • DEHP-Free, Latex-Free, and Pyrogen-Free	10 per box

