



TCBPTT063

Trifurcated Pharmacy Transfer Tubing

- + Allows simultaneous transfer from three large-volume containers, maximizing efficiency and reducing workflow interruptions.
- + Designed with precision-fit silicone tubing insert to ensure reliable performance with the Baxa Repeater® Pump for consistent and accurate fluid transfers.
- + Sterile for immediate use in pharmacy cleanrooms and sterile compounding environments.
- + Provides durability, flexibility, and optimal fluid flow control, reducing the risk of kinks and flow disruptions.



Product #	Features	Quantity
TCBPTT063	• For use with Baxa Repeater® Pump • Three Large Volume Intake Tubing • Universal Vented Spike • Silicone Tubing Insert • Male Luer Lock at Distal End • Sterile • DEHP-Free, Latex-Free, and Pyrogen-Free	10 per box

