



# TCBABFST

## Sterile Bouffant

- + Serged seams, low-linting impervious fabric, breathable polyethylene coated PET.
- + Rated for use in ISO class 5 and higher cleanroom environments.
- + Hazardous drug permeability tested.
- + One size fits all providing a snug fit.
- + Prevents exposed hair and other particles from contaminating work area.
- + Light weight and comfortable for extended wear.



| Product # | Size  | Features                                                                                                                                                                                                                                                                                                                                       | Quantity                                                          |
|-----------|-------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------|
| TCBABFST  | 21.5" | <ul style="list-style-type: none"> <li>• Elastic Opening to 21.5"</li> <li>• Breathable and Tear-Resistant</li> <li>• Latex-Free</li> <li>• Gamma Sterilized to 10<sup>-6</sup> Sterility Assurance Level</li> <li>• ASTM Permeability Testing with Hazardous Drugs and Helmke Drum Particle Testing Reports Available Upon Request</li> </ul> | Individually packed<br>10 pcs/bag<br>24 bags/case<br>240 pcs/case |

