



TCBAMKST

Sterile Face Mask

- + Excellent bacterial and particle filtration efficiency of $\geq 98\%$.
- + Fully enclosed, soft bendable nosepiece with knitted elastic ear loops.
- + External and internal layer of polypropylene spunbond, middle layer of polypropylene melt blown.
- + Recommended for ISO class 5 and higher cleanroom environments.
- + Cleanroom sterile face mask with ear loops.



Product #	Size	Features	Quantity
TCBAMKST	6.8" x 3.5"	• White 3-Ply Mask with Ear Loops • Comfortable and Breathable • Meets ASTM F21000-20 Level 3 Requirements • BFE and PFE $\geq 98\%$ • Differential Pressure < 6.0 • Splash Resistance Pressure > 160 mmHg • Flammability 16 CFR 1610: Class 1 • Latex-Free • EO sterilized to 10^{-6} Sterility Assurance Level	Individually packed 20 pcs/bag 20 bags/case 400 pcs/case

