




TCBINF6480

IV Administration Pump Set Compatible with FloGard® and Sigma Spectrum® Pumps

- + User-Friendly design simplifies the priming and loading process.
- + Made from high-quality materials to ensure durability and reliability.
- + Streamlined setup and reliable operation reduce time and effort for healthcare professionals.
- + For use in hospitals, home care, chemotherapy, hydration therapy and antibiotic administration.



Product #	Length	Drops/mL	Priming Vol.	Features	Quantity
TCBINF6480	102"	10	20mL (approx.)	<ul style="list-style-type: none"> • Universal Vented/Non-Vented Spike • PVC Fluid Path • Sterile to 10⁻⁶ SAL • 0.2µm GVS Speedflow® Filter • One Needle-Free Y-Site • Roller Clamp • Slide Clamp • Pyrogen-Free and Latex-Free • DEHP-Free PVC 	50 per box

