



TCBPTT331

Low Priming Volume Tubing

- + Minimizes drug waste, making it ideal for high-cost medications, ensuring maximum utilization and cost savings.
- + Features a precision-fit silicone tubing Insert, ensuring smooth operation and reliable performance with the Baxa Repeater® Pump.
- + Designed to meet the high-accuracy requirements of sterile compounding, reducing errors and improving workflow efficiency.
- + Supports precise medication preparation, reducing drug loss while maintaining sterility and accuracy.



Product #	Features	Qty/Unit
TCBPTT331	• Designed for High Cost Drugs • Universal Male Luer Lock • Silicone Tubing Insert • For use with Baxa Repeater® Pump • Sterile • DEHP-Free, Latex-Free, and Pyrogen-Free	10 per box

