




TCBEXT003

15" Mini-Bore Extension Set

- + Compact, highly versatile solution designed to extend the reach of primary IV tubing with minimal flow disruption.
- + Mini-bore diameter for precise fluid control, this extension set offers accuracy, durability, and compatibility across various clinical settings.
- + Features male and female luer lock connectors, compatible with standard IV setups for a leak-free and secure connection.



Product #	Length	Priming Vol.	Features	Quantity
TCBEXT003	15"	0.35mL (approx.)	• Male and Female Luer • Slide Clamp • DEHP-Free • Latex-Free • Pyrogen-Free	50 per case

