




TCBEXT004

30" Mini-Bore Extension Set

- + Provides an extended reach for primary IV tubing.
- + Enhances flexibility in patient positioning and allows for easier access to the IV site, minimizing patient discomfort and reducing the need for line adjustments during treatment.
- + Universal luer lock connectors on each end ensures a leak-free connection with standard IV lines and devices minimizing the risk of disconnection or leakage.



Product #	Length	Priming Vol.	Features	Quantity
TCBEXT004	30"	0.60mL (approx.)	• Male and Female Luer • Slide Clamp • DEHP-Free PVC • Latex-Free • Pyrogen-Free	50 per box

