




TCBEXTNFBM-7

7" Mini-Bore Bonded Needle Free Extension Set

- + Used for enhanced safety, efficiency, and reliability.
- + Lower chance of accidental disconnects or leakage.
- + Eliminates a connection point, reducing the potential for bacterial entry and contamination.
- + Eliminates the need for manual assembly, saving time in critical care settings.
- + Use for precision drug deliveries as in chemotherapy, vasoactive drugs, and sedation.
- + Pediatrics and neonates – low priming volume prevents fluid overload.



Product #	Length	Priming Vol.	Features	Quantity
TCBEXTNFBM-7	7"	0.4mL (approx.)	• Slide Clamp • DEHP-Free • Latex-Free • 7 Inch Mini-Bore • Bonded Needle-Free Valve • Rotating Male Luer Lock • Sterile	50 per box

