




TCBINF033G

IV Administration Set With GVS EasyDrop® Regulator

- + 15µm filter reduces risk of contamination, leading to safer patient outcomes.
- + GVS Flow Controller ensures IV fluids are delivered at a consistent rate, crucial for patient safety.
- + Y-site injection port allows for easy and quick administration.
- + Set offers safety, convenience, efficiency and cost effectiveness making it a valuable tool in clinical settings.



Product #	Length	Drops/mL	Features	Quantity
TCBINF033G	92"	20	• One Y-injection Site (6" From Distal End) • Male Swivel Luer Lock Connector • Universal Vented/Non-Vented Spike • Roller Clamp • 15µm Filter Located in Drip Chamber • GVS EasyDrop® Flow Regulator • DEHP and Latex free • Sterile	40 per box

