




TCBINF6419

IV Administration Pump Set with One Distal Needle-Free Y-Site

- + Compatible with FloGard® and Sigma Spectrum® Pumps.
- + User-Friendly design simplifies the priming and loading process.
- + Made from high-quality materials to ensure durability and reliability.
- + Streamlined setup and reliable operation reduce time and effort for healthcare professionals.
- + For use in hospitals, home care, chemotherapy, hydration therapy and antibiotic administration.



Product #	Length	Drops/mL	Priming Vol.	Features	Quantity
TCBINF6419	102"	10	18mL (approx.)	• Universal Vented/Non-Vented Spike • PVC Fluid Path • Sterile to 10⁻⁶ SAL • One Distal Needle-Free Y-Site • Roller Clamp • Slide Clamp • Pyrogen-Free and Latex-Free • DEHP-Free PVC	50 per box

